

Basic pH Benchtop Meters



- Simple to operate
- Automatic Temperature Compensation
- Automatic calibration
- Calibration expiration reminder
- ORP range

HI 2211 and HI 2210 are benchtop pH/ORP and °C meters. HI 2211 adds measurement for ion concentration (ISE) and Oxidation Reduction Potential (ORP) in the mV range.

These instruments also feature a reading stability indicator used during calibration and a measurement memory/memory recall function.

pH measurements for both instruments are compensated for the temperature effect manually or automatically with the HI 7662 temperature probe. These instruments are equipped with an easy-to-read LCD which shows both the primary reading and °C.

ORDERING INFORMATION

HI 2210-01 (115V), HI 2210-02 (230V), HI 2210-03 (AUS plug), HI 2211-01 (115V) and HI 2211-02 (230V) are supplied with HI 1131B pH electrode, HI 7662 temperature probe, HI 76404N electrode holder, HI 70004 pH 4.01 buffer solution sachet, HI 70007 pH 7.01 buffer solution sachet, HI 7071S electrolyte solution (30 mL), HI 700661 cleaning solution sachet, 12 VDC adapter and instructions.

ELECTRODES

HI 1043B	Use: strong acids and bases. glass body, double junction, refillable
HI 1053B	Use: emulsions. glass-body, triple ceramic junction, refillable
HI 1083B	Use: biotechnology. glass-body, open junction, refillable
HI 1131B	Use: general purpose. glass-body, ceramic junction, refillable
HI 1332B	Use: general purpose. PEI body, double junction, refillable
HI 7662	Stainless steel temperature probe with 1 m (3.3') cable

SOLUTIONS

HI 5004L	pH 4.01 buffer solution, 500 mL
HI 5007L	pH 7.01 buffer solution, 500 mL
HI 5010L	pH 10.01 buffer solution, 500 mL
HI 7061L	Electrode cleaning solution, 500 mL
HI 70300L	Electrode storage solution, 500 mL

ACCESSORIES

HI 8427	pH/mV electrode simulator
HI 931001	pH/mV electrode simulator with display
HI 76404N	Electrode holder

SPECIFICATIONS	HI 2210	HI 2211	
Range	pH	-2.00 to 16.00 pH	
	mV	-	±399.9 mV ; ±2000 mV
	Temperature	-20 to 120.0°C	
Resolution	pH	0.01 pH	
	mV	-	0.1 mV; 1 mV
	Temperature	0.1°C	
Accuracy (@20°C)	pH	±0.01 pH	
	mV	-	±0.2 mV (±399.9 mV); ±1 mV (±2000 mV)
	Temperature	±0.4 °C (excluding probe error)	
pH Calibration	automatic, one or two points with five memorized buffer values (pH 4.01, 6.86, 7.01, 9.18, 10.01)		
Temperature Compensation	automatic (with HI 7662 probe) or manual from -20.0 to 120.0°C		
pH Electrode	HI 1131B glass body pH electrode with BNC connector and 1 m (3.3') cable (included)		
Temperature Probe	HI 7662 stainless steel temperature probe with 1 m (3.3') cable (included)		
Input Impedance	10 ¹² Ohm		
Power Supply	12 VDC adapter (included)		
Environment	0 to 50°C (32 to 122°F); RH max 95%		
Dimensions	235 x 222 x 109 mm (9.2 x 8.7 x 4.3")		
Weight	1.3 kg (2.9 lbs.)		

For a complete list of Solutions and Electrodes, see the end of this section.